

Bed Head Panel

LUMIPAN



FEATURES

- a) The construction of the wall section is made of Extruded aluminium
- b) The extruded aluminium structural lower channel will include a vertical integrate equipment mounting track to support optical track mounted equipment.
- c) Optional equipments like I.V. Hook, Infusion pole, utility basket, MONitor tray etc. can be install into the track at any location along the exturded aluminium channel.
- d) Gas outlets/Switch-Socket/Monitor Socket/Nurse Call System can be fitted on middle portion of the Bed console unit after cutting the aluminium sheet as per hospital requirements.

ROBUSTRUNK



CRITRUNK



Vertical Panel

FEATURES

- a) The construction of the wall section is made of Extruded aluminium
- b) The extruded aluminium structural lower channel will include a vertical integrate equipment mounting track to support optical track mounted equipment.
- c) Optional equipments like I.V. Hook, Infusion pole, utility basket, MONitor tray etc. can be install into the track at any location along the exturded aluminium channel.
- d) Gas outlets/Switch-Socket/Monitor Socket/Nurse Call System can be fitted on middle portion of the Bed console unit after cutting the aluminium sheet as per hospital requirements.

