

BS OUTLETS

(MEDICAL GAS TERMINAL UNITS)



BRITISH STANDARD (AGSS)



BRITISH STANDARD (PANEL)



BRITISH STANDARD (WALL)

Outlets are manufactured to precision for dispensing medical gases and vacuum services reliably and safely to the secondary devices from the distribution network.

Carefully engineered, these British Standard Outlets with safe keying configurations comply with the provisions of all major standards followed internationally.

The flow has been designed to cover all the areas of a health care facility. One compatible Probe is supplied with each Outlet.

AGSS / WAGD Outlet Literature on request

Technical Information

Gas Services	Pressure	Flow Lpm
Oxygen, Nitrous Oxide, CO ₂ , Medical Air (AR 4) & Mixed Gases	58 psi (4 bar)	130
Air for Surgical Tools (AR 7) / Air Motor	102 psi (7 bar)	350
Medical Vacuum (below atmosphere)	450 mmHg (8.7 psig)	85

Ordering Information-Gas Services

BRITISH STANDARD



Safety: Gas specific Diameter Indexing (Quick Connect Type)						
Mounting	Gas Services					
	O ₂	N ₂ O	Vacuum	AR 4	AR 7	CO ₂
Wall	6013110	6013410	6013210	6013310	6013510	6013711
Panel	6013120	6013420	6013220	6013320	6013520	6013721
Bed-head	6013120	6013420	6013220	6013320	6013520	6013721
Pendant	6013140	6013440	6013240	6013340	6013540	6013741
Ceiling	6013140	6013440	6013240	6013340	6013540	6013741
Column	6013140	6013440	6013240	6013340	6013540	6013741

DISCLAIMER: Aktiv Outlets meet and exceed the requirements of relevant standards at the time of manufacture. However piping source, capacity, sizing and restrictions may prevent from attaining the stipulated value.



AKTIV TECHNOLOGIES

Sohna Flyover, Commercial Area, Opp: Devilal Park, Sohna More, Ballabhgarh,
Faridabad: 121004 (Haryana) Phone : + 91 1294062304, + 91 9811224348
E-Mail: aktivtechnologies@gmail.com, www.aktivtechnologies.com



Consult Instruction for use